

TURN YOUR VARIABLE PIPETTE INTO A FIXED PIPETTE...

The Xacto pipettor has a unique feature that decreases downtime, eases volume switching and decreases ongoing costs.

This feature allows the volume control at the top of the pipette to be switched to a fixed knob anywhere within the volume range of the pipette. When switching volumes, it's as fast as removing one knob and replacing it with another. No need to switch pipettes or change volume on the variable.

Our research shows that 76% of labs use pipettes at 2 or 3 separate volumes on a day to day basis. We've developed our adaptable volume controllers to allow the switch between these volumes faster, more efficient and more cost effective.



IN SCIENCE
TOGETHER

CODE	PRICE	DESCRIPTION	RANGE	COLOUR	DESCRIPTION
153-326	£21	5ul fixed volume controller knob	5ul	153-320	Xacto™ Pipettor, Variable Volume, 0.5 - 10ul - ISG®
153-327		10ul fixed volume controller knob	10ul		
153-328		10ul fixed volume controller knob	10ul		
153-329	£21	20ul fixed volume controller knob	20ul	153-321	Xacto™ Pipettor, Variable Volume, 2 - 20ul - ISG®
153-333	£21	30ul fixed volume controller knob	30ul	153-322	Xacto™ Pipettor, Variable Volume, 10 - 100ul - ISG®
153-334		50ul fixed volume controller knob	50ul		
153-335		100ul fixed volume controller knob	100ul		
153-336	£21	100ul fixed volume controller knob	100ul	153-319	Xacto™ Pipettor, Variable Volume, 20 - 200ul - ISG®
153-337		200ul fixed volume controller knob	200ul		
153-338	£21	500ul fixed volume controller knob	500ul	153-323	Xacto™ Pipettor, Variable Volume, 100 - 1000ul - ISG®
153-339		1000ul fixed volume controller knob	1000ul		

CONTACT OUR GLOBAL PARTNERS TODAY...



ISG® ENGLAND
Bexwell Business Park
Karoo Close
Norfolk
PE38 9GA
UNITED KINGDOM



ISG® AUSTRALIA
Po Box 1680
120 Fussell Street
Ballarat
Victoria, BT3350
AUSTRALIA



ISG® USA
ISG Westlab NA.
#4 - 17530 56A Avenue,
Surrey BC, V3S 1G3,
604 576 2776.
UNITED STATES



EXCLUSIVE PRODUCT FROM ISG®

Xacto™

PIPETTE BY ISG®

The new Xacto™ Pipettor from ISG® is the next step in pipetting revolution. Its high quality, robust European construction with a metal shaft and tip assures long lasting accuracy and limited maintenance.


This fully autoclavable pipettor features a double 'O' ring shaft for universal tip compatibility and a superb air tight seal. The Xacto™ pipette has a detachable volume mechanism for easy visibility for both right and left handlers.

The detachable volume control also allows the use of unique fixed volume knobs that can be fitted within the pipette's range rather than having to purchase a new pipettor. This allows for minimal downtime, cost effective replacement and easy recalibration.

KEY FEATURES OF THE RANGE:



THE FULL XACTO™ PIPETTOR RANGE:



CODE	DESCRIPTION	RANGE	PRICE
153-320	Xacto™ Pipettor, Variable Volume, 0.5 - 10ul - ISG®	10ul	£135
153-321	Xacto™ Pipettor, Variable Volume, 2 - 20ul - ISG®	20ul	£135
153-322	Xacto™ Pipettor, Variable Volume, 10 - 100ul - ISG®	100ul	£135
153-319	Xacto™ Pipettor, Variable Volume, 20 - 200ul - ISG®	200ul	£135
153-324	Xacto™ Pipettor, 8 Channel, Variable Volume, 20 - 200ul - ISG®	200ul	£425
153-323	Xacto™ Pipettor, Variable Volume, 100 - 1000ul - ISG®	1000ul	£135